

BLACKROCK & MICROMED

Cervello™

BLACKROCK
MICROSYSTEMS

Cervello 256

- FDA Cleared /CE Approved
- Compatible with Traditional Surface electrodes, ECoG grids, micro arrays, and all penetrating electrodes
- Network Integration
- High Definition Cameras



Cervello 256

- Up to 256 Channels
- 16 optional analog inputs channels
- Sampling rates 256 Hz to 30 kHz
- 3 separate data Streams
- System captures, processes, and analyzes (in real time) action potentials



Cervello 256

 **BLACKROCK[®]
NEUROMED**

&

**micromed**



Together the most powerful EEG system on the Market today

 **BLACKROCK
MICROSYSTEMS**

E-mail: info@blackrockmicro.com - www.blackrockmicro.com